

SAMARIT

safe patient moving & handling



Lightweight and comfortable, wherever you need it,
in the clinic, home or away ...

GLIDEBOARD

SEATED TRANSFER



CURVED GLIDEBOARD
THE GREAT
THE BIGGEST
SELFGlide
TETRAglide
ACCESSOIRES



SAMARIT

Samarit transfer systems

Transfer from bed into the wheelchair or WC, wheelchair into the motorcar etc.



THE BIGGEST

Size: 100 x 34 cm / 39 x 13 inch 2.9 kg
Order No. 476.0102

UDI/GTIN 7640144460921

This board simplifies the transfer to the patient by the nursing personnel from the centre of the bed into the wheel-chair or WC/chair (even motorcar). It is manufactured of absolutely unbreakable material (HDPE), is acid resistant, easy to clean and to disinfect. A specially developed anti-slide adhesive (blue) ensures the positioning on the bed, wheelchair, etc.



THE GREAT

Size: 80 x 33 cm / 32 x 13 inch 2.6 kg
Order No. 476.0080

UDI/GTIN 7640144460938

This board simplifies the transfer from the edge of the bed into the wheelchair or WC, wheelchair into the motorcar etc. Its composition is the same as for "The Biggest". Sit-On Trousers (accessory) is used during the transfer without clothing, e.g. onto the WC-chair. It protects the skin of the patient and is fixed with Velcro. Two solid hand grips allow the nurse to carry out the transfer in a secure and controlled manner.



CURVED GLIDEBOARD

Size: 68 x 23 cm / 27 x 9 inch 1.2 kg
Order No. 476.0012

UDI/GTIN 7640144460945

This board serves for transfer bed-wheelchair-chair etc. We developed the shape together with therapists and several tens of thousands have been sold. The shape is Design Patent Protected. The material used corresponds to the other Glideboards (HDPE).



SELFGlide

Size: 50 x 21 cm / 20 x 8 inch 0.8 kg
Order No. 476.0050

UDI/GTIN 7640144460952

The Selfglide-Glideboard is primarily used by patients who can carry out the transfer themselves due to the strength in their arms. It can also be used in the transfer with the support of the Mini-care cushion by a nurse.



TETRAGLIDE

Size: 63 x 23 cm / 25 x 9 inch 1.1 kg
Order No. 476.0863

UDI/GTIN 7640144460969

The Tetraglide Glideboard was developed specially for tetraplegics but is also increasingly popular because of its new "Fishtail". It is characterised by the ease with which the patient can be sat upon, by the manner in which it can be inserted from both sides and is utilised for transfer bed-chair-WC, etc. The material corresponds to that of the other Glideboards.

Gap < 30 cm / 12 inch, SWL* = 150 kg
Gap < 20 cm / 8 inch, SWL* = 180 kg



*SWL = Safe Weight Limit

Accessories

GLIDEBOARD

Samarit transfer systems

- Unbreakable
- Sliding surfaces
- Non slide sticker
- Lightweight
- Left/right usable
- Instruction for use
- CE mark

High quality standards
SAMARIT transfer aids have:

- Design patents
- 2 years warranty



Sit-On Monofil Sliding-Towel

Order No. 476.0200
UDI/GTIN 7640144460976

Sit-On Sliding-Towel permits easy sliding without clothes on to the Glideboards. Usually used by the patients themselves with the Curved, Selfglide and Tetraglide.

Nylon for use to glide on Glideboards to bath, shower or toilet.



Glide Belt

Order No. 450.2100	Small
UDI/GTIN 7640144461041	(72-90 cm / 28-35 inch)
Order No. 450.2200	Medium
UDI/GTIN 7640144461058	(88-106 cm / 35-42 inch)
Order No. 450.2300	Large
UDI/GTIN 7640144461065	(104-122 cm / 41-48 inch)

Helps to position the boards under heavy patients and protects the skin.

www.samarit.com

Download your language from our website or ask your local distributor.



Samarit Medical AG

Gewerbestr. 12

CH-8132 Egg, Schweiz

Tel.: +41 44 918 10 11 Fax: +41 44 918 24 39

e-mail: info@samarit.com - www.samarit.com